

# UVLI-BAGS™

PROVEN PROTECTION FOR LIGHT SENSITIVE DRUGS

Over 250 drugs require protection from ultra-violet light to ensure optimal potency. Medipak's UVLI-BAGS™ have provided that protection for over 35 years.

## ✓EASY

Lighter color green UVLI-BAGS™ make it easier to read labels and check medication levels through the bag

## ✓CONVENIENT

Our exclusive Slit-Top UVLI-BAGS™ eliminates the need to poke holes when hanging the IV

## ✓SAFE


Color code different types of drugs with green and amber bags

## ✓EFFECTIVE

Medipak manufactures the highest quality plastic film and bags. They are fully tested to meet and exceed United States Pharmacopoeia Standards

## ✓FAST SERVICE

All sizes, styles, and colors are in stock for immediate shipment\*



## EASY ACCESS

Convenient perforated tops allow for easy access

## STACKABLE

Unique three-ply cardboard design allows boxes to be stacked without crushing

## REORDER INFORMATION

Conveniently located on the front and side of every box



**FREE Samples**

**CALL NOW**

Bags are available for immediate shipment\*

**800.336.9848**

**540.667.0233**

**Fax: 540.722.1074**

9am - 5pm Eastern Time

**Questions on UV Light Protection?**

Check our FAQ's at [www.Medipak.com](http://www.Medipak.com)

For a more professional look, choose Medipak's exclusive Slit-Top UVLI-Bags for just a penny more per bag!

SLIT-TOP UVLI-BAGS™				STOCK NUMBERS			
	SIZE (WIDTH X LENGTH)	QUANTITY PER BOX	PRICE PER BOX	AMBER	GREEN	BLACK*	DARK BROWN*
For piggyback IVs (250 ml)	12,7 cm x 21,6 cm	500 Bags	\$37.00	541-SUPB-D	741-SGPB-D	841-SBPB-D	941-SDPB-D
For 1/2-liter IV bags (500 ml)	15,2 cm x 25,4 cm	1000 Bags	\$92.00	551-SUHL	751-SGHL	851-SBHL	951-SDHL
For 1-liter IV bags (1000 ml)	15,2 cm x 35,6 cm	1000 Bags	\$118.00	556-SUFL	756-SGFL	856-SBFL	956-SDFL
For 1-liter IV bags (1000 ml)	20,3 cm x 35,6 cm	500 Bags	\$71.00	561-SUXL-D	761-SGXL-D	861-SBXL-D	961-SDXL-D
For 1-liter IV bags (1000 ml)	20,3 cm x 35,6 cm	1000 Bags	\$129.00	566-SUXL-M	766-SGXL-M	866-SBXL-M	966-SDXL-M
For 3-liter IV bags (3000 ml)	25,4 cm x 45,7 cm	250 Bags	\$65.00	581-SU3L	781-SG3L	881-SB3L	981-SD3L
For 4-liter IV bags (4000 ml)	30,5 cm x 45,7 cm	250 Bags	\$73.00	586-SU4L	786-SG4L	886-SB4L	986-SD4L

REGULAR UVLI-BAGS™				STOCK NUMBERS			
	SIZE (WIDTH X LENGTH)	QUANTITY PER BOX	PRICE PER BOX	AMBER	GREEN	BLACK*	DARK BROWN*
For tablets	6,3 cm x 5,1 cm	1000 Bags	\$18.00	510-US	N/A	N/A	N/A
For large ampules	6,3 cm x 12,7 cm	1000 Bags	\$26.50	515-UM	715-GM	N/A	N/A
For syringes	6,3 cm x 21,6 cm	1000 Bags	\$37.00	520-UL	720-GL	820-BL	920-DL
For large syringes	7,6 cm x 35,6 cm	1000 Bags	\$70.00	530-ULS	730-GLS	830-BLS	930-DLS
For piggyback IVs (250 ml)	12,7 cm x 21,6 cm	500 Bags	\$32.00	540-UPB-D	740-GPB-D	840-BPB-D	940-DPB-D
For 1/2-liter IV bags (500 ml)	15,2 cm x 25,4 cm	1000 Bags	\$82.00	550-UHL	750-GHL	850-BHL	950-DHL
For 1-liter IV bags (1000 ml)	15,2 cm x 35,6 cm	1000 Bags	\$108.00	555-UFL	755-GFL	855-BFL	955-DFL
For 1-liter IV bags (1000 ml)	20,3 cm x 35,6 cm	500 Bags	\$66.00	560-UXL-D	760-GXL-D	860-BXL-D	960-DXL-D
For 1-liter IV bags (1000 ml)	20,3 cm x 35,6 cm	1000 Bags	\$119.00	565-UXL-M	765-GXL-M	865-BXL-M	965-DXL-M
For 3-liter IV bags (3000 ml)	25,4 cm x 45,7 cm	250 Bags	\$60.00	580-U3L	780-G3L	880-B3L	980-D3L
For 4-liter IV bags (4000 ml)	30,5 cm x 45,7 cm	250 Bags	\$68.00	585-U4L	785-G4L	885-B4L	985-D4L
6,3 cm Tubing	6,3 cm x 91,4 m	1 Roll	\$24.00	590-UT300	790-GT300	890-BT300	990-DT300
6,3 cm Tubing	6,3 cm x 304,8 m	1 Roll	\$50.00	590-UT1000	790-GT1000	890-BT1000	990-DT1000
6,3 cm Tubing	6,3 cm x 762 m	1 Roll	\$115.00	590-UT	790-GT	890-BT	990-DT
IV Tubing Cover 6,3 cm x 1,8 m Package of 10 precut pieces			\$9.60	590-UT6	790-GT6	890-BT6	990-DT6
IV Tubing Cover 6,3 cm x 1,8 m Box of 100 precut pieces			\$85.00	590-UT6BX	790-GT6BX	890-BT6BX	990-DT6BX

## UVLI-ZIP BAGS

SIZE (WIDTH X LENGTH)	PRICE PER 1000 (USD)	AMBER STOCK NUMBER	GREEN STOCK NUMBER	PRICE PER 100 (USD)	AMBER STOCK NUMBER	GREEN STOCK NUMBER
7,6 cm x 12,7 cm	\$84.00	AUZ35	GUZ35	\$14.00	AUZ35-C	GUZ35-C
7,6 cm x 22,8 cm	\$112.00	AUZ39	GUZ39	\$19.00	AUZ39-C	GUZ39-C
10,1 cm x 15,2 cm	\$95.00	AUZ46	GUZ46	\$16.00	AUZ46-C	GUZ46-C
12,7 cm x 20,3 cm	\$125.00	AUZ58	GUZ58	\$21.00	AUZ58-C	GUZ58-C
15,2 cm x 22,8 cm	\$135.00	AUZ69	GUZ69	\$23.00	AUZ69-C	GUZ69-C
20,3 cm x 35,6 cm	\$260.00	AUZ814	GUZ814	\$43.00	AUZ814-C	GUZ814-C
22,8 cm x 22,8 cm	\$185.00	AUZ99	GUZ99	\$31.00	AUZ99-C	GUZ99-C
22,8 cm x 30,5 cm	\$250.00	AUZ912	GUZ912	\$42.00	AUZ912-C	GUZ912-C

## CLEAR POLYBAGS (No UV Protection)

	SIZE	QUANTITY PER BOX	PRICE PER BOX	STOCK NUMBER
For tablets	6,3 cm x 5,1 cm	1000 Bags	\$12.60	610-PS
For large ampules	6,3 cm x 12,7 cm	1000 Bags	\$16.00	615-PM
For small syringes	6,3 cm x 21,6 cm	1000 Bags	\$24.00	620-PL
For medium syringes	6,3 cm x 26,6 cm	1000 Bags	\$27.60	SOP25105
For large syringes	7,6 cm x 38,1 cm	1000 Bags	\$41.00	SOP315

EMAIL ORDERS TO: [orders@medipak.com](mailto:orders@medipak.com)

FAX ORDERS TO: 540-722-1074

\*Black and Dark Brown bags & tubing covers may be subject to 1-8 week turnaround

Do you need a custom size of our UVLI-Bags? Call us for a quote.

**800.336.9848**

**540.667.0233**

**Fax: 540.722.1074**

**www.Medipak.com**

